



Size, Weight	With Single Battery: (L) 7.19" x (W) 3.26" x (H) 8.15" 4.74 lbs. / 2.15 kg. With Double Battery: (L) 7.19" x (W) 3.26" x (H) 9.03" 5.73 lbs. / 2.60 kg.
Oxygen Flow	Per setting (ml/min): 1 = 210; 2 = 420; 3 = 630; 4 = 840; 5 = 1,050; 6 = 1,260
Power	AC Power Supply: 100-240 V, 50-60 Hz (auto-sensing to allow worldwide use), DC Power for mobile use in car
Battery	Single Battery: Up to 6.5 hours Recharge Time: Up to 3 hours with AC or DC power Double Battery: Up to 13 hours Recharge Time: Up to 6 hours with AC or DC power
Noise	38 Decibels (at flow setting 2)
Warranty	5 year standard warranty*
Operation	Simple control functions and easy-to-read LCD display
Use	Designed for 24/7 use

* See OxyGo Warranty Statement for details



Stay Connected with the My OxyGo App!



Check your battery life, settings and more from your smartphone!

www.oxygo.life



Contact Your Provider for More Information:

Connect With Us



www.oxygo.life

MADE IN THE USA

Battery life statistics based on average values may vary up to +/- 10% depending on patient use.

Data accurate as of 6/19.

PN B122-0819

You Need Oxygen...
NOT TANKS.

Up To
13
HOUR
BATTERY!*



4.7 lbs. 6 Flow Settings



Bluetooth



www.oxygo.life Keep Going.™

OxyGO NEXT™

Quietest on the market
at only 38 decibels! **

7.19" x 3.26" x 8.15" 4.7 lbs.

The Next Generation of Portable Oxygen

Its ground-breaking size and battery life have never been seen before at this high of an oxygen dose, and it is setting the standard for portable oxygen concentrators. OxyGo NEXT has user-replaceable batteries, sieve beds and an easy to read LCD display with simple controls. And it can be charged at home or on the go—giving you the freedom to embrace what's NEXT!

OxyGo NEXT Comes With:

- Standard, 8-cell battery
- Nasal cannula
- Bluetooth technology
- DC car charger, AC wall charger
- Carry bag with padded, adjustable strap
- With Bluetooth technology and the **My OxyGo App**† you can troubleshoot, check battery life, oxygen purity and more with the touch of a button!



Bluetooth



My OxyGo App

Up To
13
Hour
Battery!*

6

Flow
Settings

4.7

Pounds

5

Year
Warranty



Designed to Follow You Anywhere!

OxyGo portable oxygen concentrators (POCs) are lightweight medical devices that can make the medical oxygen your doctor prescribed from the air around you. With OxyGo, you never need to worry about running out of oxygen and you don't need to carry heavy tanks!

OxyGo is different from other POC brands because the units show you exactly how much time you have left on the battery. Plus, you can change the battery easily and use it at the same time—either in the car or at home!

All OxyGo POCs are clinically validated to be used 24/7, so they can supply the medical oxygen you need at home or on the go!

* Double cell battery, setting 1. ** Flow setting 2.

† Available on OxyGo FIT POCs produced after 2/19 and OxyGo NEXT.



Single Battery
Up to 6.5 Hours!

1400-3010-8



Double Battery
Up to 13 Hours!

1400-3010-16



Carry Bag
1170-3410



Backpack
1170-3420



Slim Backpack
1170-3441



Desktop Battery Charger
1400-3030

AC Power Supply
1400-3040

DC Power Supply
1400-1050

ACCESSORIES

www.oxygo.life